


Rivertown Surgery Center
822 Farrar Drive
Conway, South Carolina
(843) 347-9587

Catheter Removal Instructions



To remove Catheter:

- Wash Hands. Place a towel under the catheter port site.
- Twist the empty syringe onto the port, as in above picture. Your catheter may or may not have a red tape.
- Allow the syringe to fill (should fill automatically). If the syringe fills completely, you may have to untwist it, expel the water, and re-twist the syringe onto the same port. Pull gently on the syringe plunger to assure that all the water is out of the catheter balloon.
- Once there is no more water filling into the syringe, gently pull catheter from the body.